

## What CelluSmooth™ offers your practice:

### Speed & Profitability

CelluSmooth is the new, innovative solution to a highly sought after procedure. Its high treatment speed produces procedure fees up to \$5,000 per hour.

### Versatility

- Offer CelluSmooth in combination with laser-assisted lipolysis. Many patients need both.
- Increase profitability by promoting combination treatments, which require little additional time.

### Superior Outcomes

- High patient satisfaction
- Minimal after-care consultation

Learn more about this breakthrough technology at [www.ScitonCelluSmooth.com](http://www.ScitonCelluSmooth.com)

## Exceptional Results

BEFORE



AFTER 6 MONTHS



Photos Courtesy of Marc J. Salzman, MD, FACS

BEFORE



AFTER 6 MONTHS



Photos Courtesy of Marc J. Salzman, MD, FACS

BEFORE



AFTER 5 MONTHS



Photos Courtesy of Hannah Vargas, MD, FACS

BEFORE



AFTER 6 MONTHS



Photos Courtesy of Marc J. Salzman, MD, FACS

## What CelluSmooth offers your patients:

### A Solution

Finally, there is an ideal solution for uneven and dimpled skin. CelluSmooth is ideal for treating these conditions and offers your patients long-lasting results for smooth looking skin.

### Get Back to Life

CelluSmooth is a minimally invasive procedure so your patients can return to everyday activities almost immediately.

### Better Results Earlier


CelluSmooth delivers results in as early as one month.

### Typically, with a single, in-office treatment

CelluSmooth offers exceptional results for a condition which effects more than 85% of women.


## CelluSmooth™ M<sup>3</sup> Technology

Three Proprietary Modes. One Unique Treatment.




### Burst Mode

A high energy pulse ablates tough, fibrous bands beneath the dermis which create dimples and valleys.



### Melt Mode

Moderate energy pulses target fat underlying bulging areas of the dermis.



### Tighten\* Mode

A gentle, low energy pulse protects the deep vascular plexus while stimulating the dermis to tighten\*.

## Superior Engineering

### Expandable Platform

Add other modules to the expandable Joule platform to grow your practice. Choose from eleven other modules offering overall skin rejuvenation, ranging from erasing mild to deep wrinkles, removing benign pigmented and vascular lesions, removing unwanted hair, firming skin and much more.



### No Consumables

We never employ life-limiting electronics or unnecessary proprietary fibers. We simply do not believe they are the right thing to do.



## A Breakthrough in Body Contouring



# What is CelluSmooth™?

CelluSmooth is a new, minimally invasive laser treatment targeting lumpy and dimpled skin.

Why CelluSmooth is unique:

### The Right Wavelength

CelluSmooth utilizes the 1319 nm wavelength which is optimal to treat:

- Dimples and valleys created by fibrous bands beneath the dermis
- Bulging areas formed by excess fat
- Thin and uneven dermis

With a penetration depth of 2-3 millimeters, combined with high power, 1319 nm is the ideal choice.

### The Power of M<sup>3</sup> Technology

M<sup>3</sup> is a proprietary technology from Sciton. It allows the delivery of 1319 nm in three distinct and unique modes with seamless transition.

- Burst
- Melt
- Tighten\*

M<sup>3</sup>'s proprietary delivery is the key to the success of the treatment.

\*Through tissue coagulation



Learn more at [www.ScitonCelluSmooth.com](http://www.ScitonCelluSmooth.com)

925 Commercial Street, Palo Alto, California 94303 • (888) 646-6999 • Sciton.com •

©2013 Sciton, Inc. All rights reserved. Sciton is a registered trademark and CelluSmooth is a trademark of Sciton, Inc. Printed in USA. 2600-029-14 Rev B

\*CelluSmooth is an application on the modular Joule system which is cleared by the FDA for the surgical incision, excision, vaporization, ablation and coagulation of soft tissue.

CelluSmooth is available for the treatment of cellulite outside of the US.

