

YOU DESERVE QUALITY...


THE **JOULE[®]X**
PLATFORM

 SCITON[®]
BECAUSE RESULTS MATTER



CUSTOM-BUILT

To Meet Your Practice Needs



1997

THE JOULE® X STORY

Since 1997, Sciton has proudly designed and manufactured in Palo Alto, CA not only the JOULE X platform, but the thousands of premium parts and components which go into it. This is just one reason for its outstanding craftsmanship, quality and longevity. JOULE X is made to order in the U.S.A. powered with proven and innovative technology advancing aesthetic and medical options for clinicians and patients worldwide.



CUSTOMIZABLE

CUSTOMIZABLE MODULAR LASER AND LIGHT PLATFORM DESIGN

Each JOULE X is custom-built with best-in-class application modules empowering your practice to fit your needs now and into the future. This value ecosystem defines upgradeability and expandability at anytime as your business grows, or as new, sensational technology becomes available like BBL® HERO.

- Broadest range of Laser & Light applications - especially popular power brand combinations
- Earn revenue fast while you ramp up - JOULE X practices' income in 6 months is typically 3X average monthly lease payment
- Sustainable approach to your business model on planned or unplanned growth
- Intelligent system health diagnostics and treatment usage monitoring with Sciton iQ

DISCOVER ROI

YOUR PRACTICE DISCOVERS IMPRESSIVE VERSATILITY YOUR PATIENTS DISCOVER THEIR BEST SKIN

JOULE X will expand your reach and attract patient populations of all skin types, ethnicities, gender, and age providing effective and desired treatment outcomes that distinguish your practice from others. Discover the powerful portfolio of applications that improve the appearance of problems or aging skin.

- High satisfaction treatment options for every patient
- Exclusive, award winning, consumer favorites only found on JOULE X
- Limitless possibilities for your practice potential

THE PERFECT SOLUTION

For Your Practice

The JOULE® X portfolio of brands is the perfect solution to expand your practice and attract new patients. Address their conditions with innovative and advanced applications.



BE
THE
PRACTICE
PATIENTS
SEEK
OUT

SKIN QUALITY

HALO®, NANOLASERPEEL™, MICROLASERPEEL®, CONTOUR TRL™, PROFRACTIONAL™, CLEARSILK®, BBL® FAMILY OF BRANDS, BBL® HERO

Textural, Wrinkles, Acne, Aging Skin

Scar Revision

PIGMENTATION

BBL® FAMILY OF BRANDS, BBL® HERO, HALO®, CLEARSILK®

Discoloration, Dullness, Uneven Tone, Age Spots

VASCULAR

BBL®, BBL® HERO, CLEARSILK®, CLEARV®

Diffuse Redness, Rosacea, Lesions

Vessels and Veins You Can and Cannot See

Unwanted Hair, Pseudofolliculitis

Onychomycosis, Plantar Warts

HAIR REMOVAL & SPECIALTIES

BBL®, CLEARTOE™

WOMEN'S HEALTH

PROFRACTIONAL™

Lesions, Scar Revision

CREATE YOUR PLATFORM

Choose best-in-class modules for best patient satisfying results and offerings, greater return-on-investment, and sustainable profitable growth. Add modules at any time.

JOULE® X Module		BBL® BroadBand Light 420-1400 - Smart Filters	Forever Bare™ 590 - 695 Smart Filters	Forever Clear™ 420, 560, 590 & 800 ST Smart Filters	Forever Young™ 560, 640, 695, 590 ST, 695 ST, 800 ST Smart Filters	SkinTyte™ 590 ST, 695 ST, & 800 ST Smart Filters	BBL® HERO - Body 515 - 1400 Smart Filters	BBL® HERO - Forever Young+™ 515 - 1400 Smart Filters	HALO® 1470 Diode & 2940 Er:YAG	Contour TRL™ 2940 Er:YAG	ProFractional™ 2940 Er:YAG	ClearSuite - ClearSilk® 1064 Nd:YAG	ClearSuite - ClearY® 1064 Nd:YAG
SKIN QUALITY	Acne	•		•									
	Acne Scars	•							•	•	•	•	
	Scar Revision								•	•	•		
	Skin Resurfacing								•	•	•		
	Skin Maintenance	•		•	•	•	•	•	•	•	•	•	•
	Skin Texture Improvement	•			•	•	•	•	•	•	•	•	•
	Wrinkles								•	•	•	•	•
PIGMENTATION	Age Spots	•			•		•	•	•	•	•	•	•
	Discoloration	•			•		•	•	•	•	•	•	•
	Dullness	•			•		•	•	•	•	•	•	•
	Uneven Tone	•			•		•	•	•	•	•	•	•
VASCULAR	Diffuse Redness - Flushing	•										•	
	Reticular Veins												•
	Rosacea	•										•	
	Spider Veins	•											•
	Telangiectasias	•			•			•				•	•
	Vascular Malformation Lesions	•											•
HAIR & SPECIALTIES	Permanent Hair Reduction	•	•										
	Pseudofolliculitis	•	•										
	Onychomycosis											•	
	Warts												•
WOMEN'S HEALTH	Lesions										•		
	Scar Revision										•		

The Sciton JOULE X Platform has multiple FDA cleared indications for its system and modules. For a complete list, contact your local representative or Sciton at info@sciton.com.

SKIN QUALITY | PIGMENTATION


Textural, Wrinkles, Scar Revision, Acne, Aging Skin
Discoloration, Dullness, Uneven Tone, Age Spots



HALO® (Hybrid Fractional Laser) — SKIN TYPE I-VI

HALO is the award-winning tunable hybrid fractional laser that delivers a non-ablative and an ablative wavelength to the same microscopic treatment zone. Combining these two wavelengths creates a strong synergistic effect, resulting in an impressive corrective solution for the return of glowing skin.

- Application driven software - treatments for *face, neck, chest, hands, body*
- LED speedometer lights ensures rapid and even treatment
- Dynamic Thermal Optimization (DTO) provides precise temperature control
- High revenue and profit generating procedure



Contour TRL



Single Spot Ablative Handpiece

CONTOUR TRL™ w/ NanoLaserPeel™, MicroLaserPeel®, and Deep Resurfacing — SKIN TYPE I-VI

Advanced resurfacing that allows you to fine tune your treatments and patient downtime. From superficial peels to deep full field resurfacing, with or without thermal coagulation, this unique precision Erbium laser scanner enables exceptional patient-tuned outcomes. Addressing a broad range of patient skin needs, from subtle to the dramatic, with superior uniformity and customizable thermal effects.

- Easy user interface which allows tunability of depth (4 µm to 200 µm), density adjustment, and coagulation thermal zone
- Deliver ablative energy in various size linear and square scanning options from a single spot to a large 30 X 30 mm area
- Comes with two single spot handpieces (2 mm and 4 mm) for targeting small areas




PROFRACTIONAL™ Therapy — SKIN TYPE I-VI


Achieve dramatic results with minimal downtime for your patients using ProFractional Therapy. It delivers ablative fractional skin resurfacing from light skin maintenance programs to more advanced treatments for exceptional scar revision.

- Easy user interface which allows tunability of depth (25 µm to 1500 µm), density adjustments, and coagulation thermal zone
- Deliver ablative energy in various size scanning options featuring sequential linear, offset and random non-sequential
- Easy visibility and maneuverability using distance stand-off

Samples of ProFractional Array Patterns



Sequential offset linear



Random non-sequential placed channels



SKIN QUALITY | PIGMENTATION

Textural, Wrinkles, Scar Revision, Acne, Aging Skin
Discoloration, Dullness, Uneven Tone, Age Spots

CLEARSilk® — SKIN TYPE I-VI



ClearSilk delivers non-ablative Nd:YAG 1064 nm wavelength to address skin's minor imperfections such as fine lines, wrinkles, acne scars and the appearance of large pores for improved luminosity, skin revitalization and subtle to refined results for all skin types, in all seasons, with zero downtime.

- Non-ablative, non-contact, comfortable
- Real-time temperature sensor for ensuring even heating of tissue
- Guided distance alignment beams for easy and fast treatments
- Audio and visual feedback

BBL® HERO and the BBL® Family of Brands — Forever Change the Way You Treat Skin



Now with the next generation BBL HERO - High Energy Rapid Output, you can quickly deliver results anywhere on the body.

The versatile and elegant designed handpiece features:
HERO - 4X Speed - 3X Peak Power - 2X Cooling

- Quickly treat large areas with Forever Body with lasting results-oriented outcomes
- Comfortable, fast, effective, low downtime
- Optimal combination and versatility of treatments
- High patient satisfaction and in-demand treatments will keep patients coming back for more!
- Clinical excellence and ultimate practice success

FOREVER CLEAR™ BBL® Phototherapy — SKIN TYPE I-V



Forever Clear BBL Phototherapy protocol features 3 specially designed target selecting Smart Filters™ to treat acne causing bacteria. They complement the full range of the comprehensive BBL phototherapy system. The CoolComfort™ precision cooling allows for a comfortable patient experience.

SKINTYTE™ BBL® Phototherapy — SKIN TYPE I-VI

SkinTyte technology features 3 specially designed infrared Smart Filters to penetrate deep into tissue. They complement the full range of the comprehensive BBL phototherapy system and are often used together. To improve the appearance of aging skin that has lost firmness, use the SkinTyte modality in a multi-pass, motion technique.



BBL® HERO - The Next Generation Broadband Light

Forever change way you treat skin from head to toe. The advanced protocols for treating aging and pigmented skin everywhere on the body produces unparalleled results with 4X Speed - 3X Peak Power - 2X Cooling.



VASCULAR

Diffuse Redness, Rosacea, Lesions, Vessels & Veins You Can and Cannot See

BBL® Vascular Phototherapy — SKIN TYPE I-V



Selected Broadband Light Smart Filters effectively and efficiently treats superficial vascularity on the face and body. Application driven software for rosacea, cherry angiomas, facial bruise, facial vessels, body vascular. Superb results in fewer treatments when compared to other IPL devices.

CLEARSilk® — SKIN TYPE I-VI



Not only does ClearSilk deliver a refined treatment for skin's imperfections, its Nd:YAG 1064 nm wavelength targets vascular chromophores. Practitioners love its ease of use and effectiveness in treating diffuse redness and rosacea.


- ClearSilk + BBL Vascular combination provides a boost to clearing microvessels
- Excellent alternative to BBL Vascular in darker skin types

CLEARV® High Power Nd:YAG — SKIN TYPE I-VI




ClearV Nd:YAG 1064 nm laser is a breakthrough solution to more effectively and comfortably treat tiny spider veins to resistant reticular vessels and vascular malformations. It features advanced engineering solution, Subsurface Vein Illumination™ (SVI™) which optimizes visualization of subsurface and feeder veins to enhance treatment. Combine with BBL for treating varying vascular pathologies.

- Easily trace individual fast-filling dynamic vessels with 3 light color and 3 light intensity options.
- Precision integrated cooling from 5 - 25 °C for safety and comfort
- Choose from 2 - 6 mm spot sizes by adjusting directly on the handpiece
- Easy user experience via intuitive application driven software



Subsurface Vein Illumination



ClearV Indication Driven User Interface

HAIR REMOVAL & SPECIALTIES

Unwanted Hair, Pseudofolliculitis, Onychomycosis, Plantar Warts



FOREVER BARE™ BBL® — SKIN TYPE I-V

A revolutionary technology in hair reduction. By sending multiple lower fluence pulses at a high repetition rate, Forever Bare BBL allows for extremely even heating of hair follicles. This minimizes missed areas that occur with traditional hair reduction devices. Use in motion or static options.

- Large volume ideal for body treatments, especially backs and legs
- Large crystal, rapid repetition rate using an in-motion technique, and precision cooling make this fast and comfortable
- Easy to use guided user interface



CLEARTOE™ — SKIN TYPE I-VI

ClearToe (utilizing the ClearSilk handpiece) delivers non-ablative Nd:YAG 1064 nm wavelength to temporarily increase a clear nail in patients with onychomycosis and for the thermal destruction of the blood supply of plantar warts.

- Non-ablative, non-contact, comfortable
- Real-time temperature sensor for ensuring even heating of target
- Guided distance alignment beams for easy and fast treatments
- Audio and visual feedback

WOMEN'S HEALTH

Lesions, Scar Revision



PROFRACTIONAL™ — SKIN TYPE I-VI

The ProFractional scanner is ideal for fractionally resurfacing external tissue for treatment of lesions and scar revision.

- Easy user interface for tunability of depth (25 μm to 1,500 μm), density adjustments, independent coagulation thermal zone adjustment
- Deliver ablative energy in various size scanning options featuring sequential linear, offset and random non-sequential
- Easy visibility and maneuverability using distance stand-off

RAVE

"If you are only going to buy one laser, you need a Sciton." **Paul Nassif, MD**

"In my office, they (Sciton) are the most reliable devices I own." **Elizabeth Tanzi, MD**

"Truly the smartest minds in all lasers." **Jill Waibel, MD**

"The JOULE Platform is a great unifier of dyschromia." **Asif Hussein, MD**

"In the loudest way...Thank you for this relationship." **Julius Few, MD**

THE SCITON WAY | Expectations and Results


SCITON CONSISTENTLY EXCEEDS EXPECTATIONS

On average Sciton devices last over a decade, in fact, most of the systems built are still in service today. Inventive, exceptional, and enduring engineering that goes beyond the commonplace. The Sciton platform has become the physician's choice for reliable high quality lasers and light based devices.



SCITON TECHNOLOGY DELIVERS RESULTS PATIENTS DESIRE

Our reputation was built on innovating technology to live up to our tagline promise of "BECAUSE RESULTS MATTER". And like Sciton, your patients' results matter, along with the long-lasting partnership you have with your patients matter. We want to build that same partnership with you.



THE PARTNERSHIP You Can Count On

At Sciton we are committed to your success. Our experienced teams have developed programs and tools for every stage of your purchasing and ownership journey.

Smart Start

5 steps to operating a thriving aesthetic practice, from launch ready, preparing your staff and office to developing your marketing plan.

Success Builder

The comprehensive Marketing, Clinical and Business success strategy offerings made up of Sciton and experienced 3rd party professionals to support your practice growth.

ScitonPro.com

You as a customer have access to this portal for essentials in building your brand with engaging and educational materials such as ready-to-use or customizable Marketing Brochures, Social Media Posts, Before and After Photos, Clinical Protocols, Treatment Videos, Training Materials, and much more!



Concierge Marketing brings industry experience, competitive knowledge and cutting edge advertising strategy to your practice. Protect your brand, be up to date with advanced insight to aesthetic treatment options. Choose from Basic and Elite Membership offerings.



Our People

Our Sciton Family is made up of YOU, our partners, our area managers, and all our employees. Our area managers are one of your most valuable resources dedicated to support your practice growth.



Service Network

Practices are busy - each time a customer calls Service or an on-site visit is scheduled our professional technicians quickly familiarize themselves with resource planning to protect your investment. Our goal is to find solutions fast to keep your platform at peak performance.



Clinical Education

Ongoing educational training is paramount in our practice support programs. We cover topics from advanced tips and techniques for impressive outcomes, to unique strategies that maximize profitability and innovative applications to expand your practice and attract new patients.



Sciton Practice Locator

Patients seek out practices who have the Sciton power brands every day. We focus on marketing campaigns with social media and key influencers to optimize the registered practice directory and reach consumers.

BECAUSE RESULTS MATTER TO YOUR PATIENTS

Visit Our Complete Gallery



HALO® | 2 months post 2 tx | courtesy of Sherrif Ibrahim, MD, PhD



CLEARV® | 3 months post 2 tx | courtesy of Peter Castillo, MD



CONTOUR TRL™ | post 1 tx | courtesy of J. David Holcomb, MD



BBL® HERO: BODY | 2 months post 1tx | courtesy of Jason Pozner, MD



BBL® HERO: BODY | post 1 tx | courtesy of Patrick Bitter, Jr., MD, FAAD



PROFRACTIONAL™ | post 4 tx | courtesy of Jeff Walding, MD

BECAUSE RESULTS MATTER TO YOUR PATIENTS

Visit Our Complete Gallery



CLEARSilk® | 1 month post 4 tx | courtesy of Rebecca Gelber, MD



CLEARSilk® + FOREVER YOUNG™ BBL® | 6 years of yearly routine tx total of 8 tx | courtesy of Antonio Campo Voegeli, MD



CLEARV® + BBL® | 2 tx | courtesy of Antonio Campo Voegeli, MD



HALO® + BBL® | 4 weeks post 1 tx | courtesy of Dermatology Associates Plymouth Meeting



CONTOUR TRL™ + PROFRAC TIONAL™ + SKINTYTE™ | post 1 tx | courtesy of Michael Lin, MD



BBL® + MICROLASERPEEL® | 2 weeks post 1 tx | courtesy of Jacque P. LeBeau, MD

SCITON OWNS SKIN...

OWN THE SCITON EXPERIENCE WITH JOULE® X AND mJOULE™

Designing and building results-oriented technology for treating skin is our passion. Now Sciton offers two expandable platforms that have exceptional applications. Both platforms feature award winning BBL HERO.

HALO
CONTOUR TRL
PROFRACTIONAL
NANOLASERPEEL
MICROLASERPEEL
CLEARSIK
CLEARV
BBL FAMILY OF BRANDS
BBL HERO FACE & BODY

Legacy
JOULE X is the high performance, multi-module platform that continues to be the preferred choice for the most selection of laser and BBL applications. The workstation that quickly becomes the workhorse of practices. It is your perfect path for expansion of services.



BBL FAMILY OF BRANDS
BBL HERO FACE & BODY
SKINTYTE
MOXI

mJOULE, an affordable luxury platform, offers the same great innovative technology as JOULE X. Debuting on the mJOULE is MOXI™, a non-ablative fractional laser together with next generation BBL, BBL HERO. These two impressive technologies offer solutions that appeal to the younger prejuvenation patient, those new to lasers or who are looking to improve their skin care regimen.

YOU DESERVE
BETTER
YOU DESERVE
SCITON



Corporate Headquarters

925 Commercial Street, Palo Alto, California 94303
Phone: (888) 646-6999 • Email: info@sciton.com
www.sciton.com

Sciton Worldwide

United Kingdom, Australia, Canada, Japan

©2021 Sciton, Inc. All rights reserved. Sciton, JOULE, MicroLaserPeel, BBL, HALO, ClearSilk, and ClearV are registered trademarks and Contour TRL, NanoLaserPeel, ProFractional, ClearToe, Forever Bare, Forever Clear, Forever Young, BBL HERO Body, SkinTyte, MOXI, and mJOULE are trademarks of Sciton, Inc.
Printed in USA

2600-029-04 Rev. O