THE SCITON EDGE

SCITONPRO.COM

Complimentary Sciton pro website:

- Digital practice marketing assets
- Available 24/7

SMART START™ GUIDE

5 steps to operating a thriving practice.

SUCCESS BUILDER™

Access to experienced third-party pros and tools to support your practice growth.

- Clinical education
- Marketing tools
- Business development

CONCIERGE MARKETING

Bring industry experience, competitve knowledge and cutting edge advertising strategy to your practice.

- Customized social media, email & newsletters
- Personalized support for website, B&As and menu of services
- Hands-on help through Zoom calls & presence at in-person events.







SPECIFICATIONS: CLEAR SUITE 1064 nm

Wavelength	1064 nm Nd:YAG
Fluence	Up to 400 J/cm ²
Pulse Width	Up to 200 ms
Repetition Rate	Up to 20 Hz
Spot Size	ClearSilk: 7 mm ClearV: 2-6 mm
Skin Cooling	ClearV: Adjustable, 5 °C to 25 °C

SPECIFICATIONS: BBL® HERO™

Wavelength	420 – 1400 nm
Filters	420 nm, 515 nm, 560 nm, 590 nm, 640 nm, 695 nm, 800 nm
Output	Up to 30 J/cm ²
Pulse Width	Up to 200 ms
Repetition Rate	Up to 4 pulse-per-second
Cooling Method	Continuous thermoelectric sapphire plates
Skin Cooling	Adjustable from 0 - 30 °C in 1 °C increments
Spot Size	15 x 45 mm
Finesse Adapters	15 x 15 mm square, 11 mm & 7 mm round







©2021 Sciton, Inc. All rights reserved. Sciton, JOULE, HALO, BBL, ClearV, ClearSilk, BroadBand Light (BBL), and Forever Young are registered trademarks and Finesse Adapters, Smart Filters, HERO, Resurfacing, MOXI and SVI are trademarks of Sciton, Inc.









Clear Suite THE FAMILY OF BRANDS

Clear Clear \$1/k.





Clear Suite

Sciton is excited to provide you with a comprehensive vascular solution for your practice combining laser and light to treat a myriad of indications all year long on all skin types.

This not only further elevates the JOULE® platform as the most customizable laser and light platform in the industry, but also ensures a continuous revenue stream and great Return on **Investment for YOU!**

> I have used 1064 Nd:YAG for years and I find the ClearSilk® in combination with BBL® HERO™ to be a superior choice for treating rosacea. The combination boosts the effect on the control of flushing and inflammatory flare-ups. No other system offers me the versatility that my JOULE provides.

> > Antonio Voegeli-Campo, MD, PhD

Sciton's ClearV® is the new gold standard for vascular treatments! It's unique Subsurface Vein Illumination allows me to clearly see veins below the surface and even treat feeder veins I can't see for a more complete and precise treatment. Coupled with integrated cooling and adjustable spot sizes, ClearV makes the treatment much more comfortable for the patient and easier for the provider.

Antonio Voegeli-Campo, MD, PhD

As a dermatologist I appreciate having an entire arsenal of modalities in one platform to treat multiple skin conditions. The Clear Suite in combination with BroadBand Light offers me an unmatched capability to treat vascular lesions and improve the signs of sun damaged skin.

Sherrif F. Ibrahim, MD

The Clear Suite 1064 Nd:YAG applications - ClearV and ClearSilk - combined with BroadBand Light (BBL) offer you the unique capability to treat vascular lesions and revitalize skin from a single platform.

APPLICATIONS		ClearV 1064 nm	ClearSilk 1064 nm	BBL HERO
SKIN QUALITY	Acne			•
	Acne Scars			•
	Pigmented Lesions		•	•
	Scar Revision		•	
	Skin Maintenance		•	•
	Wrinkles		•	
	Fine lines		•	•
VASCULAR	Vascular Lesions	•	•	•
	Telangiectasia	•		•
	Reticular/Spider Veins	•		
	Diffuse Redness		•	
HAIR & SPECIALTIES	Permanent Hair Reduction			•
	Onychomycosis (ClearToe)		•	
	Warts	•		
ОТНЕВ	Skin Types	All	AII	I-V
	Anesthetic (numbing) Required	None	None	None
	Days of Downtime	None	None	None
	All Seasons	•	•	•

CLEAR SUITE: The Family of Brands | 2 $\textbf{1} \mid \textbf{CLEAR SUITE: The Family of Brands}$



Treat the Veins You Can and Can't See!

ClearV® is the preferred choice for vascular treatments. It effectively and safely treats visible vessels, broken capillaries, spider veins, vascular lesions, and skin discoloration. Gently treats undesirable, visible micro-vessels and resistant, unwanted subsurface vessels by applying targeted heat to stimulate the body's natural processes to clear away.

Technology

VascuZoom

Enables quick and easy spot size changes directly on the handpiece. This enables the user to switch effortlessly between 2-6 mm spot sizes.

Subsurface Vein Illumination (SVI™)

SVI technology is designed to give the provider enhanced visibility to see subsurface and feeder veins to treat more thoroughly and effectively. Push button controls on the handpiece offer the ability to change the lighting color, level of brightness, and aiming beam intensity.





Clear

Indication Driven Software

Allows to easily determine treatment range and safe start parameter settings.

Adjustable Integrated Cooling

Sapphire tip provides optimal cooling pre, during and post as well as bulk cooling for patient comfort. The adjustable cooling range of 5-25 °C is ideal for treating a variety of vascular conditions safely, efficaciously, and efficiently.

Indications



REDNESS, i.e. DIFFUSE REDNESS AND ROSACEA



FACIAL VESSELS



PORT WINE STAINS





CHERRY & SPIDER ANGIOMAS

VENOUS

LAKES

ClearV Benefits

Better visualization and control through SVI™

All Fitzpatrick skin types

No disposables

No patient downtime – treat all year long

Easily combine with BBL HERO and other Sciton modalities





ClearV | 3 months post 2 tx | courtesy of courtesy of Peter Castillo, MD





ClearV | 8 weeks post 1 tx | courtesy of courtesy of Jason Pozner, MD





ClearV | 4 weeks post 2 tx | courtesy of Antonio Voegeli-Campo, MD, PhD

3 | CLEAR SUITE: The Family of Brands | 4



Boost Skin Revitalization

The ClearSilk® treatment works by deposition of heat, stimulating fibroblast activity and the body's natural repair process resulting in a decrease in fine lines and wrinkles. ClearSilk is the treatment of choice for all skin types. Patients have no downtime and may resume daily activities immediately following treatment.

Technology

Built-in Temperature Sensing

Easily visualize the skin temperature as you are treating with the built-in sensing gauge that indicates temperature ranges between 32-48 °C.



Audible Feedback

The ClearSilk handpiece and JOULE® system deliver audible alarms to provide real time feedback on the number of pulses delivered.



Visual Distance Guide

Enable user to maintain consistent treatment distance to target for uniform heating to ensure reliable results.

Indications



REDNESS, i.e. DIFFUSE REDNESS AND ROSACEA



FINE LINES AND WRINKLES



PIGMENTED LESION



ACNE SCARS

ClearSilk Benefits

All Fitzpatrick skin types

Easily delegatable

Zero downtime

No anesthetic or disposables required

Comfortable with minimal patient discomfort

Non-ablative, and non-contact

4-6 treatments recommended to achieve optimal results





ClearSilk | 1 month post 4 tx | courtesy of Rebecca Gelber, MD





ClearSilk | 6 months post 3 tx | courtesy of Yonick Plastic Surgery





ClearSilk | 6 months post 3 tx | courtesy of Peter Castillo, MD

5 | CLEAR SUITE: The Family of Brands | 6



The next generation BroadBand Light® from Sciton, BBL HERO is finally here. BBL now with High Energy Rapid Output™, or HERO, can quickly deliver results anywhere on the body.

- 4x Speed 3x Peak Power 2x Cooling
- Quickly treat large areas with BBL HERO Body with outstanding outcomes
- Comfortable, fast, effective, low downtime
- Optimal combination and versatility of treatments
- High patient satisfaction and in-demand treatments will keep patients coming back!
- Clinical excellence and ultimate practice success.

BBL HERO | 2 tx | courtesy of Jason Pozner, MD

BBL HERO Benefits

Treat broad range of skin types, I-V, and conditions

Combine with ClearV®, ClearSilk®, Resurfacing[™], MOXI[™] and HALO[®] to offer best in class treatments

Fast, effective, versatile, and cost effective

Predictable and consistent

User friendly

Technology

Now with CoolComfort™ Technology

The most powerful cooling capacity allowing the user to set and maintain cooling from 1-30 °C in 1 °C increments for fast, comfortable and uniform treatments.



Uniquely designed with each lamp supplying half the energy during a pulse resulting in an extended lifetime of the BBL handpiece.

Indications

MILD TO MODERATE **ACNE VULGARIS** 420

PIGMENTED LESIONS, INCL. LENTIGOS (AGE SPOTS, SOLAR LENTIGOS (SUN SPOTS)

VASCULAR LESIONS, SUCH AS DIFFUSE REDNESS, ROSACEA. AND TELANGLECTASIAS. VASCULAR COMPONENT OF SCAR REVISION

PORT WINE STAINS, HEMANGIOMA, VASCULAR 590 MALFORMATIONS, HAIR REMOVAL, DARK SKIN FILTERS

HAIR REMOVAL 640 DARK SKIN FILTERS

HAIR REMOVAL DARK SKIN FILTERS

590ST SKIN FIRMING

695ST SKIN FIRMING

800ST SKIN FIRMING



BBL HERO | post 1 tx | courtesy of Jason Pozner, MD





BBL HERO | 1 month post 3 tx | courtesy of Steven Swengel, MD





BBL HERO | 3 tx | courtesy of of Patrick Bitter, Jr., MD, FAAD

7 | CLEAR SUITE: The Family of Brands

COMBINATION TREATMENT & SPECIALTIES

ClearV® and ClearSilk® are highly versatile and can easily be combined not only

with BBL®, but with many other Sciton modalities.

Clear / + FBBL

ClearV and BBL together provide a great tool

and valuable solution for complex vascular conditions, such as Port Wine Stains.





ClearV + BBL | 2 tx | courtesy of Antonio Campo Voegeli, MD, PhD





ClearV + BBL | 6 months post 8 tx | courtesy of John Cooksey, MD

Clear Silk and BBL together are a powerful





ClearSilk + BBL | 1 month post 2 tx | courtesy of Antonio Campo-Voegeli, MD, PhD

solution to treat rosacea and diffuse redness effectively and lastingly.





ClearSilk + BBL | post 2 tx | courtesy of Antonio Campo-Voegeli, MD, PhD

In addition, the ClearV and ClearSilk handpieces can be used to treat other nail and the nail bed resulting in a decrease of nail fungus and improve nail





7 months post 3 tx | courtesy of Marla Quintana, RN

conditions, such as warts and onychomycosis. The laser energy heats the growth and appearance.





8 months post 1 tx | courtesy of Mary Beth Mudd, MD