



JOULE™

Modular, Multi-Application Laser and Light Platform for
Aesthetic, Medical and Surgical Applications



joule THE WORLD'S MOST EXPANDABLE PLATFORM



BUILD THE PERFECT SYSTEM FOR YOUR PRACTICE

The Sciton JOULE™ is the world's only expandable platform to offer an unprecedented array of laser and light wavelengths through any of its three distinct delivery modes: arm, fiber and BroadBand Light™. As a result, JOULE allows practitioners to provide the widest range of aesthetic, medical and surgical procedures from a single platform.

CHOOSE FROM A COMPREHENSIVE RANGE OF TREATMENTS



SKIN | Resurfacing: Hybrid, Full Field, Fractional

HALO™ NanoLaserPeel™, MicroLaserPeel® Contour TRL™ and ProFractional™



SKIN | Phototherapy™

BBL™, Forever Young BBL™ and Forever Clear™



VASCULAR

BBL™, ClearScan YAG™, ClearScan ALX™ and Pro-V™



BODY

ALLURA™, CelluSmooth™ and SkinTyte™



HAIR & SPECIALITIES

Forever Bare BBL™, ClearScan YAG™ and ClearScan ALX™

HOW TO CREATE YOUR SYSTEM

MIX & MATCH
Combine treatment modalities to create a completely customized system.

NEED TO FILL A VOID?
The **JOULE PLATFORM** is the most flexible system to satisfy a unique applications mix.

JUST GETTING STARTED?
Create a **specialized system** that has the ability to grow with your practice.

NEED IT ALL?
Try the **fully loaded system**, which offers the most comprehensive range of treatments.

ARM DELIVERY



HALO™ (Hybrid Fractional Laser)
Halo is the world's first tunable Hybrid Fractional Laser combining ablative and non-ablative wavelengths in a single pass to maximize results and control downtime.



ProFractional Therapy – ProFractional™ and ProFractional-XC/III™
Achieve dramatic results with minimal downtime for your patients using Sciton's ProFractional Therapy for fractional skin resurfacing with optional coagulation.



Contour TRL™ with MicroLaserPeel® and NanoLaserPeel™
From superficial to deep peels, the unique Contour Tunable Resurfacing Laser enables exceptional, patient-tuned outcomes. Address a broad range of patient needs, from the subtle to the dramatic, with superior uniformity and customizable thermal effects.



ClearScan ALX™
ClearScan ALX provides superior patient comfort, effective and trusted results, and fast treatment times for permanent hair reduction, vascular and benign pigmented lesions, and wrinkles.



ClearScan YAG™
An effective solution for permanent hair reduction and the treatment of virtually all pigmented and vascular lesions, ClearScan YAG is the gold-standard technology for all skin types.



ThermaScan™
Invigorate any skin type and tanned skin with ThermaScan, a non-ablative, infrared laser that effectively and precisely treats wrinkles, acne, acne scars, and other skin blemishes.



ClearSense™ Accessory
Extend the versatility of the JOULE platform by offering real-time temperature feedback and optimum laser delivery to safely and effectively treat a variety of conditions that include skin and nail rejuvenation.

FIBER DELIVERY



ALLURA™
Enabling impressive body sculpting results, ALLURA is best-in-class for laser-assisted lipolysis. This ultra high-powered, multiple wavelengths laser melts unwanted fat and provides tissue coagulation for a tighter, slimmer body.



TempASSURE™ Accessory
Safely deliver high power laser-assisted lipolysis with accurate, internal, real-time temperature monitoring provided by the TempASSURE accessory for ALLURA.



CelluSmooth™
A breakthrough in body contouring and new minimally invasive laser treatment targeting lumpy and dimpled skin. 1319 nm is the ideal choice with a penetration depth of 2-3 millimeters combined with high power.



Pro-V™
Treat venous insufficiency with excellent precision and minimal pain, bruising, and swelling. Endovenous laser ablation with Pro-V technology provides unmatched control of vein coagulation and ensures a comfortable treatment with consistent results.

BBL DELIVERY



BBL™ with Forever Young BBL™ and Forever Clear BBL™
With a single versatile handpiece and interchangeable smart filters, BBL is a comprehensive phototherapy system and is an excellent choice to use for vascular lesions, to treat acne, permanent hair reduction, and much more. Its elegant design, unsurpassed lamp life, and breadth of applications give it superior value.



SkinTyte™
Quickly achieve the appearance of firmer skin and high patient satisfaction with SkinTyte. This economical and innovative light energy device is ideal for rapid selective heating of dermal collagen and for treating large areas.



Forever Bare BBL™
A revolutionary technology in hair removal. By sending multiple lower fluence pulses at a high repetition rate, Forever Bare BBL allows for extremely even heating of hair follicles and minimizes missed areas that occur with traditional hair removal devices.

