



BBLINK™ Goggles and BBL® HERO™ – The Perfect Pair

We've Raised The Bar For Pulsed Light Eye Protection



True Block



Increase Treatment Volume



Reduce Eyestrain
and Headaches

Proprietary shutter design automatically matches the speed of any BBL® application - HERO™, Forever Young™, Forever Bare™, and SkinTyte™. We do the BBLink-ing for you. Consistent and reliable protection from every flash, every time.

- Own the best - Stand out from the crowd
- Headaches and eye strain from basic IPL protection are a thing of the past
- Better visibility for precise placement of the BBL sapphire for perfect results every time!

"I wouldn't want to perform a BBL pulsed light treatment without them now! They're terrific!" S. Barbarich, RN


"Blinkers are a huge improvement over existing BBL glasses. I don't want to use HERO or BBL without them." Jason Pozner, MD

BBLINK™

BBLINK Goggles - Features

First of their kind, BBLINK™ Goggles feature proprietary technology which connects directly with your BBL® system and receives shutter close commands from the system computer slightly before the BBL flash occurs. This allows them to keep up with the high-speed pulsing of HERO™, while still allowing you to see everything but the flash. No other shutter goggles on the market today match the flash speeds or syncing capability. BBLINK works with all modes including HERO, Forever Young™, Forever Bare™, and SkinTyte™. Blocks every flash, every time.


- Artfully designed with integrated CE marked IPL protective goggles, BBLINK offers True Block
- Easy to use and easy on your eyes
- Quick and simple pairing with your JOULE®, mJOULE™ System*
- Multiple goggles can be synced to the same system and used during the same treatment
- Can last an average full day of treatments
- Proprietary design blinks so fast you almost wonder if the BBL is pulsing at all



Experience the perfect pair for yourself. Your eyes will THANK you!

Let the goggles do your BBLINK-ing

You Deserve Better | You Deserve Sciton

 Proudly Made in the U.S.A.